Cross Reference Chart **APILUS**



		JUNIOR 3G	SENIOR3G	xCell Pur	x Cell Pro
	13.56 MHz	✓	✓		
Frequency	27.12 MHz			✓	✓
Thermolysis	Thermolysis (13.56 MHz HF current working in tenths of a second)	✓			
	Flash (13.56 MHz HF current working in hundredths of a second)	✓	✓		
	MicroFlash (13.56 MHz HF current working in thousandths of a second)		✓		
	PicoFlash (27.12 MHz HF current working in thousandths of a second)			✓	✓
	MeloFlash (27.12 MHz HF current working in hundredths of a second)			✓	✓
	Pp/vd function (allows the application of two reduced MicroFlash or PicoFlash pulses in two different spots inside the hair follicle)		✓	✓	~
	Synchro (synchronized application of multiple ultrafast 27:12 MHz HF micropulses)			optional	✓
	MultiPlex (slow thermolysis with a MicroFlash or PicoFlash pulse at the end)		~	✓	~
Blend	Pulsing Blend (continuous GLV current with 13.56 MHz HF pulses)	✓	✓		
	PicoBlend (continuous GLV current with 27.12 MHz HF pulses)			✓	✓
	OmniBlend (combination of continuous GLV and HF currents)	✓	✓	✓	✓
	MultiBlend (combination of continuous GLV and HF currents with a MicroFlash or PicoFlash pulse at the end)		✓	✓	✓
	SynchroBlend (continuous GLV current with multiple 27.12 MHz HF micropulses)			optional	✓
	EvoluBlend (simultaneous application of a gradually increasing GLV current and 27.12 MHz HF current)			optional	✓
Electrolysis (continuous GLV current)		✓	✓	✓	✓
Printer					✓
Tolerance test		✓	✓	✓	✓
Moisture test*		✓	✓	✓	✓
Number of programs (covering each body area, hair size and probe type)		512	1080	19800	29700
Number of personalized programs			60	60	100
Cataphoresis mode (soothing treatment)		✓	✓	✓	✓
Anaphoresis mode		✓	✓	✓	✓
VitaPhase lonomassage* device		optional	optional	optional	✓
Screen	Regular (80 characters)	✓	✓		
	Touchscreen (800 × 480 pixels)			✓	✓
Key to calculate U.L. (units of lye)			✓	✓	✓
Key to calculate E L (energy level)			✓	✓	✓
Operation mode without footswitch		✓	✓	✓	✓
Management software and 1000 client files					✓
Medical records, appointment calendar and statistics					✓
Smart card system					✓
USB port and updates				✓	✓
Number of probe holder outlets		1	2	2	2
Accessory test			✓	✓	✓
Current protection and stability monitor		✓	✓	✓	✓
International power supply				✓	✓

^{*} not available in the USA